


Disrupt, homogenize,
dissociate or lyse cells,
tissues and plant matter



Stomp your samples flat

In seconds. Right in the tube.

FEATURES

- ▶ Can accommodate 5 mL and microcentrifuge tubes.
- ▶ Powerful lever action makes even the toughest samples easy to crush.
- ▶ Quick and easy. Stomping only takes a few seconds.
- ▶ Disposable pestle shield prevents cross-contamination.



Your problem samples won't know what hit them

With the Next Advance Stomper™, you can crush tough material in seconds, with no cross-contamination.

Use it with the Bullet Blender® homogenizer to process samples like soybeans, shells, and corn!

Collect



Stomp



Homogenize



Step 1: Place your sample in the tube.

Step 2: Stomp it.

Step 3: Add grinding beads and homogenize in the Bullet Blender®. Your sample will be a fine powder in as *little as 5 minutes*.

For More Information, Contact: